HOW TO ATTACH YOUR

KNEE GEL CUSHION FOR TKW KNEE WALKER

Item # TKW-GC



Follow these SIMPLE STEPS

CONTACT US!
We're here to help

info@novajoy.com

1 800 557 6682

novajoy.com

1 Pull one side of the strap down **(see Figure 1).**



(Figure 1)

2 Secure the strap by pulling it underneath the knee post **(see Figure 2).**



(Figure 2)

3 Do the same for the other strap **(see Figure 3).**



(Figure 3)

4 Ensure that both straps are securely attach before each use **(see Figure 3).**



(Figure 4)

To Clean the Gel Core

Remove the cover and spot clean the gel core or wipe down with a warm moist cloth. After cleaning allow the gel to dry completely before replacing the cover and using.

To clean the Cover

To remove the cover, fold the cushion and peel the cover back and off of the gel cushion. Machine-wash in cold water and air dry ONLY. Daily maintenance can be done easily using a moist sponge or towel.

