



novada

349

Parts Breakdown
For models: 349B, 349R
Serial # Starting "NV"



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on December 4, 2018

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your transport chair corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).


If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!



Figure 1

Main Body Components




Handle and Folding Mechanism Components




Detachable Arm Components



Footrest Assembly



Front And Rear Wheel Assembly Components





Item Number 349B, 349R

349 TRANSPORT CHAIR

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
VT-AT100	ANTI-TIPPER (PAIR) BLACK <i>SEE SN REFERENCE BELOW</i>	PAIR
PVT-33215	ARM PAD DESK LENGTH	EACH
VT-AS100	ARM SUPPORT BRACKET	EACH
PVT-WA100	AXLE FOR FRONT WHEEL	EACH
PVT-WA102	AXLE FOR REAR WHEEL	EACH
CH-31929	BEARING & SPACER FOR FRONT WHEEL	SET
PVT-WB101	BEARING FOR 8" REAR WHEEL	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
VT-33201	BLUE LEFT DETACHABLE ARM	EACH
VT-33202	BLUE RIGHT DETACHABLE ARM	EACH
VT-BL101	BRAKE 8" LEFT WHEEL	EACH
VT-FC101	COVER FOR FOOTREST	EACH
VT-BR101	BRAKE FOR 8" RIGHT WHEEL	EACH
PVT-AC102	COVER FOR REAR WHEEL	EACH
PVT-FC100	COVER FOR FRONT FORK SWIVEL SCREW	EACH
V-33004	CRADLE FOR FRONT SEAT	EACH
VT-AP100	DESK ARM PANEL	EACH
VT-FP102	FOOT PLATE FOR LEFT FOOTREST	EACH
VT-FP101	FOOT PLATE FOR RIGHT FOOTREST	EACH
VT-F100B	FOOT REST FOR BLUE	PAIR
VT-F100R	FOOT REST FOR RED	PAIR
VT-FF100	FORK FOR FRONT WHEEL	EACH
VT-SG100	GLIDE FOR SEAT	PAIR
VT-HG100	HAND GRIP	PAIR
VT-HS101	HEEL LOOP FOOT REST	EACH
V33001	KNOB FOR ADJUSTABLE FOOTREST	SET
PVT-FM100	LEFT FOLDING MECHANISM	EACH
PVT-RL100	LEVER TO RELEASE FOLD DOWN ARM	EACH
PVT-AN100	NUT FOR AXLE - FRONT/REAR	EACH
VT-P101	PIN FOR FOLDING MECH	EACH
VT-P100	PIN FOR FOOTREST	PAIR
PVT-PL101	PLUG FOR FRAME	PAIR
PVT-PL100	PLUG FOR SEAT BAR	PAIR
VT-LF106B	PUSH HANDLE LEFT - BLUE	EACH
VT-RF106B	PUSH HANDLE RIGHT- BLUE	EACH



Item Number 349B, 349R

349 TRANSPORT CHAIR


SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
VT-LF106R	PUSH HANDLE LEFT - RED	EACH
VT-RF106R	PUSH HANDLE RIGHT- RED	EACH
S5000AH	RECEIVER FOR FLIP BACK ARM-REAR	EACH
VT-33203	RED LEFT DETACHABLE ARM	EACH
VT-33204	RED RIGHT DETACHABLE ARM	EACH
PVT-RL101	RELEASE LEVER ON FLIP BACK ARM	EACH
PVT-FM101	RIGHT FOLDING MECHANISM	EACH
VT-TK101	RUBBER TOE KICK	EACH
PVT-CS102	SCREW FOR CROSSBAR (INCLUDES NUT, SPACER & WASHERS)	EACH
V33003	SCREW TO ATTACH DESK ARM PAD	PAIR
PVT-HS100	SCREW TO ATTACH HEEL LOOP TO FOOT REST	EACH
VT-33205	SCREW TO ATTACH SIDE PANEL	EACH
VT-SU101	SCREWS FOR BACK UPHOLSTERY	SET
VT-33206	SCREW FOR FRONT SEAT CRADLE	EACH
VT-SU100	SCREW FOR SEAT UPHOLSTERY	SET
CT-SB100	SEAT BELT	EACH
P40003	SILENCER FOR ANTI-TIPPERS	EACH
VT-WS101	SPACER FOR REAR WHEEL AXLE	EACH
VT-33208	SPACER FOR SWIVEL FRONT FORK	EACH
CK-31903	TIP FOR BRAKE HANDLE	EACH
VT-UB104	UPHOLSTERY 19" BACK	EACH
VT-US100	UPHOLSTERY 19" SEAT	EACH
PVT-FW100A	WHEEL 8" GREY FRONT	EACH
PVT-RW100A	WHEEL 8" GREY REAR	EACH

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>POST JULY 2018 MODELS</i>
349B	SN BEFORE: NV349BHE0103
349R	SN BEFORE: NV349RHF0001

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE