



Cruiser De-Light

Parts Breakdown
For models: 4010BL, 4010RD
Serial # Starting "NT"



CARE • QUALITY • INNOVATION



THE NOVA PROMISE

Revised on January 28, 2020

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see **Figure 1**).

If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!




Figure 1


Upper Assembly Components




Lower Assembly Components



Rear Brake Assembly Components



Front Wheel Assembly Components





Item Number 4010

Cruiser De-Light

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
V420201	AXLE & NUT FOR FRONT WHEEL	EACH
P40107	AXLE FOR REAR WHEEL	EACH
4010B	BASKET	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P40006	BEARINGS FOR 5"/6"/8" WHEEL	SET
P40010	BEARING PLASTIC FOR FRONT SWIVEL FORK	EACH
P40105LB	BLUE LEFT FRONT LEG EXT. W/ 5in. WHEEL	EACH
P40105RB	BLUE RIGHT FRONT LEG EXT. W/5" WHEEL	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
P-COVER	BRAKE COVER	EACH
P40106L	BRAKE SHOE FOR 5" LEFT WHEEL	EACH
P40106R	BRAKE SHOE FOR 5" RIGHT WHEEL	EACH
P401013	COVER FOR BRAKE SHOE	PAIR
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
P401010	CROSSBAR - BLUE	EACH
P401012	CROSSBAR - RED	EACH
P401011	FOLDING RELEASE PULL KNOB	EACH
P42022	FORK FOR 5" FRONT	EACH
P40104	FRAME CONNECTOR	EACH
B15-C	HAND BRAKE WITH STEEL CABLE	EACH
P40103B	LEFT REAR LEG EXT. W/ 5in. WHEEL - BLUE	EACH
P40103R	LEFT REAR LEG EXT. W/ 5in. WHEEL - RED	EACH
P40011	NUT FOR FORK STEM	EACH
P40016	NUT FOR FRONT/REAR WHEEL AXLE (BAG OF 5)	SET
CH-902608	PUSH PIN & SPRING (SET OF 4)	SET
P40105LR	RED LEFT FRONT LEG EXT. W/ 5in. WHEEL	EACH
P40105RR	RED RIGHT FRONT LEG EXT. W/ 5in. WHEEL	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH
P40104B	RIGHT REAR LEG EXT. W/ 5IN. WHEEL-BLUE	EACH
P40104R	RIGHT REAR LEG EXT. W/ 5in. WHEEL-RED	EACH
P40103	SCREW & NUT TO ATTACH HAND BRAKE	PAIR
P40101	SCREW & NUT FOR BRAKE SHOE	EACH
P40110	SCREW & NUT TO ATTACH SEAT TO FRAME	EACH
P40007	SCREW TO ADJUST CABLE TENSION	EACH
P40035	SCREW TO ATTACH FORK TO FRAME	EACH
P40100	SCREW WITH BARREL CLAMP FOR BRAKE CABLE	EACH



Item Number 4010

Cruiser De-Light

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
P40102	SEAT WITHOUT PAD	EACH
P42058	SEAT PAD	EACH
P40109	SPRING FOR BRAKE SHOE	EACH
P42019	WHEEL 5"	EACH

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE