



# nova

# Cruiser Deluxe

Parts Breakdown  
*For models: 4202BK, 4202BL, 4202GN, 4202PL, 4202RD  
Serial # Starting "NT"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

*Revised on November 11, 2019*

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

## An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1


# Hand Brake Mechanism




# Upper Assembly Components



## Seat & Backrest



## Folding Mechanism & Accessories



# Wheel & Lower Assembly Components





Item Number 4202

**Cruiser Deluxe**

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
P42004	ANTI-RATTLE SEAT CLAMP	PAIR
P40021	AXLE FEMALE FOR REAR WHEEL	EACH
P42043	AXLE FOR FRONT WHEEL	EACH
P40017	AXLE REAR WHEEL - MALE	EACH
P42193BK	BACKREST BLACK FLIP-UP	EACH
P42193BL	BACKREST BLUE FLIP-UP	EACH
P42193GN	BACKREST GREEN FLIP-UP	EACH
P42193PL	BACKREST PURPLE FLIP-UP	EACH
P42193RD	BACKREST RED FLIP-UP	EACH
C4203B	BASKET FOR WALKERS	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P40006	BEARINGS FOR 5"/6"/8" WHEEL (BAG OF 5)	SET
P-CABLE	BRAKE CABLE (PAIR)	PAIR
P40027	BRAKE SHOE	PAIR
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V42093	FOAM FOR FLIP-UP BACKREST	EACH
P42191FM	FOLDING MECHANISM (FITS BOTH LEFT AND RIGHT SIDES; SCREWS INCLUDED).	EACH
KP4218F	FORK ASSEMBLY FOR 8" FRONT WHEEL	EACH
P40033	FORK FOR 8" FRONT WHEEL W/STEM	EACH
P40031	FORK FOR 8" REAR WHEEL	EACH
P43003	GRIP - LEFT ANATOMICAL	EACH
P43004	GRIP - RIGHT ANATOMICAL	EACH
B14-C	HAND BRAKE WITH STEEL CABLE	EACH
P421512	HANDLE SILENCER FOR 4215/4202/	EACH
P421514	HANDLE SPACER	EACH
V420902	HARDWARE FOR FOLDING MECHANISM	EACH
KP4202	KIT FOR NARROW HUB 4202 WHEEL	EACH
H-1002	LOCK N LOAD	EACH
P40015	NUT COVER FOR WHEEL AXLES	EACH
P40028	NUT FOR AXLE, FOLDING MECH. AN A-FRAME	EACH
P40011	NUT FOR FORK STEM	EACH
P40016	NUT FOR FRONT/REAR WHEEL AXLE	EACH
P42058	SEAT PAD	EACH
P42041	PLUG FOR SLANTED FRAME FOLDING	EACH



Item Number 4202


### Cruiser Deluxe

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH
V42064	SCREW & NUT FOR BACKREST HINGE	EACH
P42007	SCREW FOR FOLDING MECH	EACH
V42009	SCREW SET FOR FRAME	EACH
P40007	SCREW TO ADJUST CABLE TENSION	EACH
P40030	SCREW TO ATTACH BRAKE SHOE TO REAR FORK	EACH
P40035	SCREW TO ATTACH FRONT FORK TO FRAME	EACH
C40006	SCREW TO ATTACH HAND BRAKE	EACH
P40029	SCREW TO ATTACH REAR FORK TO FRAME	EACH
P40001	SCREWS TO ATTACH PLASTIC SEAT	EACH
P42155	SEAT	EACH
P42042	SPACER FOR REAR WHEEL AXLE	EACH
P42191	STEEL HANDLE (LEFT)	EACH
P42192	STEEL HANDLE (RIGHT)	EACH
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT	EACH
P42003	U-SHAPE SEAT ATTACHMENT	EACH
P42009	WASHER FOR FOLDING MECHANISM	EACH
P42010	WASHER PLASTIC CONTOURED FOR FOLDING MECH	EACH
P42085	WHEEL 8"	EACH

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

CARE • QUALITY • INNOVATION



THE NOVA PROMISE