



# GetGo Classic

Parts Breakdown

*For models: 4202CBK, 4202CBL, 4202CDB, 4202CGN,  
4202CPK, 4202CPL, 4202CRD  
Serial # Starting "NV"*



CARE • QUALITY • INNOVATION



THE NOVA PROMISE

*Revised on December 20, 2018*

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

## An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see **Figure 1**).

If you need further assistance, please call our customer service team at **(800) 557-6682**.

**Thank you!**



**Figure 1**

# Hand Brake Mechanism

P-CABLE



P-COVER



P40007



C40006



P40047




B14-C



## Upper Assembly Components


V42080



V42010



V42079



V42094



P40003



H-1002




V42089




V42009




## Seat & Backrest



## Folding Mechanism & Accessories



# Wheel & Lower Assembly Components





Item Number 4202C

**GetGo Classic**

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V420201	AXLE & NUT FOR FRONT WHEEL	EACH
V42097	AXLE FOR REAR WHEEL	EACH
V42099	BACKREST BLACK FLIP-BACK	EACH
V42065	BACKREST BLUE FLIP-BACK	EACH
V42102	BACKREST FLIP SKYBLUE	EACH
V42067	BACKREST GREEN FLIP-BACK	EACH
V42098	BACKREST PINK FLIP-BACK	EACH
V42020	BACKREST PURPLE FLIP-BACK	EACH
V42066	BACKREST RED FLIP-BACK	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P40006	BEARINGS FOR WHEEL	BAG5
P42037	BLACK STRAP CROSSBAR	EACH
P-CABLE	BRAKE CABLE (PAIR)	EACH
V42087	BRAKE SHOE	PAIR
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	EACH
V42093	FOAM FOR FLIP-UP BACKREST	EACH
P42006	FOLDING CROSSBAR	EACH
P42038	FOLDING MECHANISM	EACH
KP4202F	FORK ASSEMBLY FOR FRONT WHEEL <i>PRE FEB. 2015</i>	EACH
P42065	FORK FOR FRONT WHEEL	EACH
B14-C	HAND BRAKE WITH STEEL CABLE	EACH
V42010	HAND GRIP ANATOMICAL	PAIR
V42079	HANDLE FOR LEFT SIDE	EACH
V42080	HANDLE FOR RIGHT SIDE	EACH
V420902	HARDWARE FOR FOLDING MECHANISM	EACH
H-1002	LOCK N LOAD	EACH
P40015	NUT COVER FOR WHEEL AXLES	EACH
P42012	NUT FOR BRAKE SHOE SCREW	EACH
P40028	NUT FOR FOLDING MECH	EACH
P40011	NUT FOR FORK STEM <i>PRE FEB. 2015</i>	EACH
P40016	NUT FOR FRONT/REAR WHEEL AXLE	BAG5
V42094	PLUG FOR SLANTED FOLDING FRAME	EACH
V4218N	POUCH W/INNER POCKET	EACH



Item Number 4202C

**GetGo Classic**

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH
V42064	SCREW & NUT FOR BACKREST HINGE	EACH
P42015	SCREW FOR BRAKE SHOE	EACH
P42007	SCREW FOR FOLDING MECH. & SLANTED FRAME PLUG	EACH
V42077	SCREW FOR PADDED SEAT	SET4
V42009	SCREW SET FOR FRAME	EACH
P40007	SCREW TO ADJUST CABLE TENSION	EACH
P40035	SCREW TO ATTACH FORK TO FRAME <i>PRE FEB. 2015</i>	EACH
C40006	SCREW TO ATTACH HAND BRAKE	EACH
V42007	SEAT CLAMP	EACH
V42084	SEAT- PADDED W/POCKET	EACH
P40003	SILENCER FOR ADJ. HANDLE	EACH
P42013	SPRING FOR BRAKE (PAIR)	PAIR
V42011	STOPPER FOR SEAT (PAIR)	PAIR
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT KNOB	EACH
P40012	WASHER FOR AXLE	EACH
P42009	WASHER FOR FOLDING MECHANISM CONTOURED	EACH
P40023	WASHER FOR FRONT & REAR AXLE	EACH
P42010	WASHER PLASTIC CONTOURED FOR FOLDING MECHANISM	EACH
P42016	WASHER PLASTIC CONTOURED REAR WHEEL	EACH
V42052D	WHEEL 6" GROOVED	EACH
P40014	WHITE PLASTIC WASHER FOR BRAKE	EACH

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>PRE FEB. 2015 MODELS</i>
4202CBL	<b>SN BEFORE:</b> NV4202CBLDM0221
4202CGN	<b>SN BEFORE:</b> NV4202CGNDM0246
4202CPK	<b>SN BEFORE:</b> NV4202CPKDN0131
4202CRD	<b>SN BEFORE:</b> NV4202CRDDN0141

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

CARE • QUALITY • INNOVATION



THE NOVA PROMISE