



nova

Zoom 20

Parts Breakdown
*For models: 4220BK, 4220BL, 4220DB, 4220GN, 4220PK,
4220PL, 4220RD
Serial # Starting "NV"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on February 5, 2019

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (**see Figure 1**).

If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!



Figure 1

Hand Brake Mechanism

P-CABLE



P-COVER



P40007



C40006



P40047



B14-C



Upper Assembly Components

V42080



V42010



V42079



V42094



P40003



H-1002




V42089




V42009




Seat & Backrest



Folding Mechanism & Accessories



Wheel & Lower Assembly Components





Item Number 4220

Zoom 20

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V42081	AXLE FOR FRONT WHEEL	EACH
V42082	AXLE FOR REAR WHEEL- FEMALE	EACH
P40017	AXLE REAR WHEEL - MALE	EACH
V42099	BACKREST BLACK FLIP-BACK	EACH
V42065	BACKREST BLUE FLIP-BACK	EACH
V42067	BACKREST GREEN FLIP-BACK	EACH
V42098	BACKREST PINK FLIP-BACK	EACH
V42020	BACKREST PURPLE FLIP-BACK	EACH
V42066	BACKREST RED FLIP-BACK	EACH
V42102	BACKREST SKYBLUE FLIP-BACK	EACH
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P40006	BEARINGS FOR 5"/6"/8" WHEEL	EACH
P42037	BLACK STRAP FOR FOLDING MECHANISM	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
V42087	BRAKE SHOE	EACH
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V42093	FOAM FOR FLIP-UP BACKREST	EACH
V42088	FOLDING MECHANISM	EACH
V42001	FORK FRONT	EACH
V42075	FORK FRONT <i>PRE FEB. 2015</i>	EACH
B14-C	HAND BRAKE WITH STEEL CABLE	EACH
V42010	HAND GRIP ANATOMICAL (PAIR)	PAIR
V42079	HANDLE BAR LEFT	EACH
V42080	HANDLE BAR RIGHT	EACH
H-1002	LOCK N LOAD	EACH
P40028	NUT FOR AXLE, FOLDING MECH. AND A-FRAME	EACH
P42012	NUT FOR BRAKE SHOE SCREW	EACH
P40011	NUT FOR FORK STEM <i>PRE FEB. 2015</i>	EACH
V42085	NUT FOR FRONT AXLE	EACH
V42094	PLUG FOR SLANTED FOLDING FRAME	EACH
V4218N	POUCH W/INNER POCKET	EACH
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH



Item Number 4220

Zoom 20

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V42064	SCREW & NUT FOR BACKREST	EACH
V42077	SCREW TO ATTACH PADDED SEAT	SET
P42015	SCREW FOR BRAKE SHOE	EACH
P42007	SCREW FOR FOLDING MECH. & SLANTED FRAME PLUG	EACH
V42009	SCREW SET FOR A-FRAME	SET
P40007	SCREW TO ADJUST CABLE TENSION	EACH
P40035	SCREW TO ATTACH FORK TO FRAME <i>PRE FEB. 2015</i>	EACH
C40006	SCREW TO ATTACH HAND BRAKE	EACH
V42007	SEAT CLAMP	EACH
V42084	SEAT- PADDED	EACH
P40003	SILENCER FOR ADJ. HANDLE	EACH
V42090	SPACER CONTOURED FOR REAR WHEEL	EACH
V42091	SPACER FOR REAR WHEEL AXLE	EACH
P42013	SPRING FOR BRAKE (PAIR)	PAIR
V42002	STEM & WASHERS FOR V42001 FORK	EACH
V42011	STOPPER FOR SEAT (PAIR)	PAIR
V42089	TRIANGLE SHAPED HEIGHT ADJUSTMENT	EACH
P42009	WASHER FOR FOLDING MECHANISM	EACH
P42010	CONTOURED WASHER FOR FOLDING MECHANISM	EACH
V42074D	WHEEL 8"	EACH
P40014	WHITE PLASTIC WASHER FOR BRAKE SCREW	EACH

SERIAL NUMBER (SN) REFERENCE

ITEM NUMBER	<i>PRE FEB. 2015 MODELS</i>
4220BK	SN BEFORE: NV4220BKDM0041
4220BL	SN BEFORE: NV4220BLDM0151
4220GN	SN BEFORE: NV4220GNDM0001
4220RD	SN BEFORE: NV4220RDDM0101

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE