



Traveler

Parts Breakdown
*For models: 4900BK, 4900BL, 4900DB, 4900PL, 4900RD
Serial # Starting "CH"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on November 21, 2023

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame near the rear, right wheel (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1


Hand Brake Mechanism



Upper Assembly Components



Lower Assembly Components



Fork & Wheel Assembly

CH-490005



CH-490006



FRONT

JN-490022




CH-490001



CH-490003



CH-490010



PLEASE REFERENCE
SN CHART


CH-490013



CH-490004



REAR



CH-490002



	Item Number 4900	
	Traveler	
	SERIAL # BEGINNING WITH "CH"	
ITEM NUMBER	DESCRIPTION	UNIT
CH-490001	AXLE & NUT FRONT	EACH
CH-490002	AXLE & NUT REAR	EACH
JN-490005	BAG	EACH
JN-490006	BASKET	EACH
CH-490003	BEARING & SPACER FOR WHEEL	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
JN-423621	BRAKE CABLE HOLDER	EACH
CH-490004	BRAKE SHOE (SET)	SET
CH-490005	CAP FOR SWIVEL FORK	EACH
JN-490004	COVER FOR CABLE	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
JN-490022	BEARING (PAIR) FOR FRONT SWIVEL FORK	PAIR
JN-490016	FOLDING MECHANISM	EACH
CH-490006	FORK SWIVEL & NUT	EACH
JN-490001	GRIP (PAIR)	PAIR
JN-490003	HAND BRAKE WITH CABLE	EACH
JN-490017	HANDLE BAR LEFT	EACH
JN-490018	HANDLE BAR RIGHT	EACH
JN-420439	REFLECTOR FOR HAND BRAKE	PAIR
C40006	SCREW TO ATTACH BRAKE	EACH
JN-490002	SCREW WITH KNOB (PAIR) FOR HT. ADJUSTMENT	PAIR
CH-490007	SILENCER FOR FRAME	PAIR
CH-490008	SUPPORT BRACE	EACH
CH-490009	TRAY FOR 4900 BASKET	EACH
CH-490010	WHEEL 8"	EACH
CH-490013	NEW WHEEL 8" FOR 4900	EACH

CHANGE EFFECTIVE	
ITEM	STARTING SERIAL #
4900BK	CH4900BKPK0001
4900BL	CH4900BLPK0001
4900DB	CH4900DBPK0001
4900PL	CH4900PLPK0001
4900RD	CH4900RDPK0001

CONTACT US!
We're here to help

✉ info@novajoy.com
☎ 1 800 557 6682
novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE