



# novada

## 332

Parts Breakdown  
For models: 332B, 332R  
Serial # Starting "NV"



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on December 4, 2018

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

## An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your transport chair corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).


If you need further assistance, please call our customer service team at **(800) 557-6682**.

**Thank you!**



Figure 1

## Main Body Components



## Hand Brake And Rear Components







## Detachable Arm Components



## Footrest Assembly



# Front And Rear Wheel Assembly Components





Item Number 332B, 332R

**332 TRANSPORT CHAIR**

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
VT-AT100	ANTI-TIPPER (PAIR) BLACK <i>SEE SN REFERENCE BELOW</i>	PAIR
PVT-33215	ARM PAD DESK LENGTH	EACH
VT-33209	ARM SUPPORT BRACKET LEFT	EACH
VT-33210	ARM SUPPORT BRACKET RIGHT	EACH
PVT-WA100	AXLE FOR FRONT WHEEL	EACH
PVT-WA101	AXLE FOR REAR WHEEL	EACH
CH-31929	BEARING & SPACER FOR FRONT WHEEL	SET
PVT-WB102	BEARING FOR 12" WHEEL	SET
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
VT-COVER	BLACK CABLE COVER	EACH
VT-33201	BLUE LEFT DETACHABLE ARM	EACH
VT-33202	BLUE RIGHT DETACHABLE ARM	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
VT-BL100	BRAKE FOR 12" LEFT WHEEL	EACH
VT-BR100	BRAKE FOR 12" RIGHT WHEEL	EACH
VT-FC101	COVER FOR FOOTREST	EACH
PVT-FC100	COVER FOR FRONT FORK SWIVEL SCREW	EACH
PVT-AC101	COVER FOR REAR WHEEL AXLE	EACH
V-33004	CRADLE FOR FRONT SEAT	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
VT-AP100	DESK ARM PANEL	EACH
VT-FP102	FOOT PLATE FOR LEFT FOOTREST	EACH
VT-FP101	FOOT PLATE FOR RIGHT FOOTREST	EACH
VT-F100B	FOOT REST FOR BLUE	PAIR
VT-F100R	FOOT REST FOR RED	PAIR
VT-FF100	FORK FOR FRONT WHEEL	EACH
VT-SG100	GLIDE FOR SEAT	PAIR
B14-C	HAND BRAKE WITH STEEL CABLE	EACH
V42010	HAND GRIP ANATOMICAL	PAIR
VT-HS101	HEEL LOOP FOR FOOT REST	EACH
V33001	KNOB FOR ADJUSTABLE FOOTREST	EACH
PVT-AN100	NUT FOR FRONT AXLE	EACH
PVT-AN102	NUT FOR REAR AXLE	EACH
VT-P100	PIN FOR FOOTREST	PAIR
PVT-PL101	PLUG FOR FRAME	PAIR



Item Number 332B, 332R

## 332 TRANSPORT CHAIR

SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
PVT-PL100	PLUG FOR SEAT BAR	PAIR
S5000AH	RECEIVER FOR FLIP BACK ARM-REAR	PAIR
VT-33203	RED LEFT DETACHABLE ARM	EACH
VT-33204	RED RIGHT DETACHABLE ARM	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	EACH
PVT-RL101	RELEASE LEVER ON FLIP BACK ARM	EACH
VT-TK101	RUBBER TOE KICK	EACH
PVT-CS101	SCREW FOR CROSSBAR (INCLUDES NUT, SPACER & WASHERS)	SET
PVT-BS101	SCREW FOR HAND BRAKE	PAIR
VT-SU101	SCREW FOR BACK UPHOLSTERY	SET
VT-33206	SCREW FOR FRONT SEAT CRADLE	EACH
VT-SU100	SCREW FOR SEAT UPHOLSTERY	SET
PVT-BS100	SCREW SET FOR WHEEL BRAKE	EACH
P40007	SCREW TO ADJUST CABLE TENSION	EACH
V33003	SCREW TO ATTACH DESK ARM PAD	EACH
PVT-HS100	SCREW TO ATTACH HEEL LOOP	EACH
VT-33205	SCREW TO ATTACH SIDE PANEL	EACH
CT-SB100	SEAT BELT	EACH
P40003	SILENCER FOR ADJ. HANDLE	EACH
VT-WS100	SPACER FOR REAR WHEEL AXLE	EACH
VT-33208	SPACER FOR SWIVEL FRONT FORK	EACH
CK-31903	TIP FOR BRAKE HANDLE	EACH
VT-UB102	UPHOLSTERY 22" BACK	EACH
VT-US102	UPHOLSTERY 22" SEAT	EACH
V33002A	WHEEL 12" GREY REAR	EACH
PVT-FW100A	WHEEL 8" GREY FRONT	EACH

### SERIAL NUMBER (SN) REFERENCE

ITEM NUMBER

*POST AUG. 2018 MODELS*

332B

SN AFTER: NV332BHF0001

### SERIAL NUMBER (SN) REFERENCE

ITEM NUMBER

*POST SEPT. 2018 MODELS*

332R

SN AFTER: NV352RHF0001

CONTACT US!  
We're here to help

✉ [info@novajoy.com](mailto:info@novajoy.com)  
☎ 1 800 557 6682  
[novajoy.com](http://novajoy.com)

CARE • QUALITY • INNOVATION



THE NOVA PROMISE