



nova

STAR HD

PETITE

Parts Breakdown
*For models: 4279BL, 4279PL, 4279RD
Serial # Starting "NV"*




CARE • QUALITY • INNOVATION



THE NOVA PROMISE

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

Revised on February 11, 2020

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (**see Figure 1**).

If you need further assistance, please call our customer service team at **(800) 557-6682**.

Thank you!




Figure 1

Hand Brake Mechanism




Upper Assembly Components



Frame & Seat Components



Wheel & Lower Assembly Components






Item Number 4279

STAR HD PETITE


SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V428817	AXLE FRONT	EACH
V428818	AXLE REAR	EACH
V427911	BACKREST BLUE	EACH
V427913	BACKREST PURPLE	EACH
V427910	BACKREST RED	EACH
V428825	BAG (INCLUDES PLASTIC LINING)	EACH
P40006	BEARINGS FOR WHEEL	BAG5
P-CABLE	BRAKE CABLE (PAIR)	PAIR
V428810	BRAKE SHOE (SET)	SET
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V42093	FOAM FOR FLIP-UP BACKREST	EACH
V427912	FOLDING MECHANISM INCLUDES HARDWARE	EACH
V428815	FRAME ATTACHMENT	EACH
V428816	FRAME ATTACHMENT HARDWARE	EACH
V428827	FRAME INSERT/ LOCK & LIFT HOLDER	EACH
B15-C	HAND BRAKE WITH STEEL CABLE	EACH
V42010	HAND GRIP ANATOMICAL (PAIR)	PAIR
V427907	HANDLE BAR LEFT	EACH
V427908	HANDLE BAR RIGHT	EACH
V427901	HANDLE/LEG ASSEMBLY LEFT	EACH
V427902	HANDLE/LEG ASSEMBLY RIGHT	EACH
V427810	LEG FRONT ASSEMBLY	EACH
V427903	LEG REAR	EACH
V428824	LOCK & LIFT	EACH
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	PAIR
P40103	SCREW & NUT (PAIR) TO ATTACH HANDBRAKE	PAIR
P40007	SCREW TO ADJUST CABLE TENSION	EACH
V42077	SCREWS TO ATTACH PADDED SEAT	SET
V427905	SEAT	EACH
V42007	SEAT CLAMP	EACH
V428820	SILENCER FRONT LEG FRAM	EACH
V428819	SILENCER REAR LEG FRAME	EACH
V428813	STEM & WASHERS	EACH
V42011	STOPPER FOR SEAT (PAIR)	PAIR

	Item Number 4279	
	STAR HD PETITE	
	SERIAL # BEGINNING WITH "NV"	
ITEM NUMBER	DESCRIPTION	UOM
V428809	TOE KICK (PAIR)	PAIR
V427906	WHEEL 8"	EACH

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE