



nova

STAR HD

Parts Breakdown
For models: 4289BL, 4289PL, 4289RD
Serial # Starting "NV"



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on February 28, 2024

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1

Hand Brake Mechanism




Upper Assembly Components



Frame & Seat Components



Wheel & Lower Assembly Components






Item Number 4289

STAR HD


SERIAL # BEGINNING WITH "NV"

ITEM NUMBER	DESCRIPTION	UOM
V428817	AXLE FRONT FOR 4288/4278 SN: NV	EACH
V428818	AXLE REAR (SET) FOR 4288/4278 SN: NV	EACH
V428907	BACKREST BLUE	EACH
V428909	BACKREST PURPLE	EACH
V428908	BACKREST RED	EACH
V428906	BAG FOR 4289 (INCLUDES PLASTIC LINING)	EACH
P40006	BEARINGS FOR WHEEL	BAG5
P-CABLE	BRAKE CABLE (PAIR)	PAIR
V428810	BRAKE SHOE (SET)	SET
P-COVER	CABLE COVER BLACK	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
V42093	FOAM FOR FLIP-UP BACKREST	EACH
V428901	FOLDING MECHANISM	EACH
V428815	FRAME ATTACHMENT	EACH
V428816	FRAME ATTACHMENT HARDWARE	EACH
V428827	FRAME INSERT/ LOCK & LIFT HOLDER	EACH
B15-C	HAND BRAKE WITH STEEL CABLE	EACH
P-TKW1229	HAND GRIP ANATOMICAL (PAIR)	PAIR
V428910	HANDLE BAR LEFT	EACH
V428911	HANDLE BAR RIGHT	EACH
V428902	HANDLE/LEG ASSEMBLY LEFT	EACH
V428903	HANDLE/LEG ASSEMBLY RIGHT	EACH
V428905	LEG REAR	EACH
V428824	LOCK & LIFT	EACH
V42076	RECEIVER FOR REMOVABLE BACKREST	EACH
P40047	REFLECTOR FOR BRAKE HANDLE	PAIR
P40103	SCREW & NUT (PAIR) TO ATTACH HANDBRAKE	PAIR
P40007	SCREW TO ADJUST CABLE TENSION	EACH
V42077	SCREWS TO ATTACH PADDED SEAT	SET
V427905	SEAT	EACH
V428828	SEAT CLAMP	EACH
V428820	SILENCER FRONT LEG FRAM	EACH
V428819	SILENCER REAR LEG FRAME	EACH
V428813	STEM & WASHERS	EACH
V428814	STEM & WASHERS FOR FORK	EACH
V42011	STOPPER FOR SEAT (PAIR)	PAIR

	Item Number 4289	
	STAR HD	
	SERIAL # BEGINNING WITH "NV"	
ITEM NUMBER	DESCRIPTION	UOM
V428809	TOE KICK (PAIR)	PAIR
V427906	WHEEL 8"	EACH

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

CARE • QUALITY • INNOVATION



THE NOVA PROMISE