



# nova

## 7200

Parts Breakdown  
*For models: 7200L, 7200LE, 7201L, 7201SE  
Serial # Starting "YH"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

*Revised on July 14, 2023*

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

## An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your wheelchair corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame near the front wheel (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!



Figure 1


# Main Body Components



## Rear Frame Components




## Footrest & Lower Frame Components






# Front Wheel Assembly

SEE SERIAL NUMBER REFERENCE AT THE BOTTOM



# Rear Wheel Assembly






Item Number 7200L/7200LE/7201L/7201LE

## 7200 Lightweight Wheelchair

SERIAL # BEGINNING WITH "YH"

| ITEM NUMBER | DESCRIPTION                                | UOM  |
|-------------|--|------|
| YU-710016   | ANTI-RATTLE CLAMP AND SCREW                | EACH |
| YH-716004   | ANTI-TIPPER BLACK ADJ (PAIR)               | PAIR |
| YU-506019   | ARM PAD DESK ARM WITH SCREWS               | EACH |
| YU-506052   | ARM PAD FULL LENGTH WITH SCREW             | EACH |
| YU-506015   | AXLE SET FOR 8" WHEEL                      | SET  |
| YU-506013   | BEARING & SPACER SET FOR 24" WHEEL         | SET  |
| YU-506012   | BEARING & SPACER SET FOR 8" WHEEL          | SET  |
| YU-506041   | BRAKE LEFT                                 | EACH |
| YU-506042   | BRAKE RIGHT                                | EACH |
| YU-506011   | CAP FOR FOOTREST (PAIR)                    | PAIR |
| S5100SC     | CRADLE FOR REAR SEAT                       | EACH |
| YH-716002   | DESK ARM LEFT ALUM H/T                     | EACH |
| YH-716003   | DESK ARM RIGHT ALUM H/T                    | EACH |
| W-7001      | ELEVATING LEG RESTS                        | PAIR |
| YU-710007   | FOOT PLATE (PAIR) WITH HEEL LOOP           | PAIR |
| YH-716001   | FOOTRESTS BLACK (PAIR)                     | PAIR |
| YU-506045   | FORK FOR FRONT WHEEL                       | EACH |
| YH-716008   | FULL ARM LEFT ALUM H/T                     | EACH |
| YH-716009   | FULL ARM RIGHT ALUM H/T                    | EACH |
| YU-506039   | GRIP (PAIR)                                | PAIR |
| YU-506056   | HAND RIM FOR 24" WHEEL - INCLUDES 9 SCREWS | EACH |
| YH-716005   | HANDLE BAR H/T                             | EACH |
| YU-710019   | HEEL LOOP                                  | EACH |
| YU-506060   | HUB CAP FOR 24" WHEEL                      | PAIR |
| YU-710023   | KNOB & SCREW FOR HANDLE BAR                | EACH |
| YU-710026   | LEVER TO RELEASE FLIP BACK ARM             | EACH |
| YU-506054   | PLUG (PAIR) FOR FOOT REST                  | PAIR |
| YU-710021   | RECEIVER & CRADLE FRONT LEFT               | EACH |
| YU-710022   | RECEIVER & CRADLE FRONT RIGHT              | EACH |
| YU-710024   | RECEIVER REAR FOR DESK ARM                 | EACH |
| YU-506055   | SCREW FOR ARM PAD (2 PCS)                  | PAIR |
| YU-710020   | SCREW FOR CENTER OF CROSSBAR               | EACH |
| YU-710018   | SCREW FOR HEEL LOOP                        | EACH |
| YU-506034   | SEAT BELT FOR WHEELCHAIR                   | EACH |
| YU-710017   | SPACER CONTOURED                           | PAIR |
| S5100CS2    | SUPPORT FOR CROSS BAR                      | EACH |
| S5100SB     | SWIVEL BEARING (PAIR)                      | PAIR |
| YU-506051   | TIP FOR BRAKE                              | EACH |


|  | Item Number 7200L/7200LE/7201L/7201LE |      |
|---|---------------------------------------|------|
|   | <b>7200 Lightweight Wheelchair</b>    |      |
|   | SERIAL # BEGINNING WITH "YH"          |      |
| ITEM NUMBER   | DESCRIPTION                           | UOM  |
| YU-710003   | UPHOLSTERY BACK FOR 7200 NYLON        | EACH |
| YU-710004   | UPHOLSTERY SEAT FOR 7200 NYLON        | EACH |
| YU-506005   | 8" WHEEL FRONT                        | EACH |
| YH-716006   | 24"WHEEL REAR WITH CAP                | EACH |

| SERIAL NUMBER FOR REFERENCE |                                 |
|-----------------------------|---------------------------------|
| 7200L                       | SN STARTING FROM: YH7200LNJ0413 |

| ITEM NUMBER | DESCRIPTION                       | UOM  |
|-------------|-----------------------------------|------|
| YU-506065   | AXLE SET FOR 8" WHEEL             | SET  |
| YU-506064   | BEARING & SPACER SET FOR 8" WHEEL | SET  |
| YU-506066   | FORK FOR FRONT WHEEL              | EACH |
| YU-506063   | WHEEL 8" FRONT                    | EACH |

**CONTACT US!**  
We're here to help

 [info@novajoy.com](mailto:info@novajoy.com)

 1 800 557 6682

[novajoy.com](http://novajoy.com)

CARE • QUALITY • INNOVATION



THE NOVA PROMISE