



no^vo^va

TKW-12

Parts Breakdown
*For models: TKW-12
Serial # Starting "NT"*



CARE • QUALITY • INNOVATION




THE NOVA PROMISE

Revised on October 11, 2023

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com

An important note from



Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your Knee Walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker on the side of the frame (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1



Steering And Brake Components



Front And Rear Wheel Components



Knee Pad Platform And Accessories






Item Number TKW-12

THE KNEE WALKER

SERIAL # BEGINNING WITH "NT"

ITEM NUMBER	DESCRIPTION	UOM
P-TKW1205	AXLE SET (PAIR) FOR FRONT AND REAR WHEEL	PAIR
P-TKW1215	BARREL SCREW SET (PAIR) FOR BRAKE	PAIR
TKW-BASKET	BASKET	EACH
P-TKW1405	BEARING FOR LOWER STEERING POST	PAIR
TKC-1023	BEARING FOR WHEEL	EACH
P-TKW1207	BOLT SET FOR FOLDING HINGE	EACH
P-TKW1201	BRAKE LEFT	EACH
P-TKW1227	BRAKE PAD (PAIR)	PAIR
P-TKW1202	BRAKE RIGHT	EACH
P-TKW1406	C-CLAMP STEERING	EACH
P-TKW1214	CABLE W/COVER FOR BRAKE	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
P-TKW1229	GRIP (PAIR) FOR TKW	PAIR
P-TKW1213	HAND BRAKE W/CABLE	EACH
B15-C	HAND BRAKE WITH STEEL CABLE SEE SN REFERENCE BELOW	EACH
V42010	HAND GRIP ANATOMICAL SEE SN REFERENCE BELOW	PAIR
P-TKW1225	HANDLE BAR	EACH
P-TKW1301	HANDLE BAR SEE SN REFERENCE BELOW	EACH
P-TKW1228	HANDLE BAR SEE SN REFERENCE BELOW	EACH
P-TKW1224	HEIGHT ADJUSTMENT LEVER (ORANGE)	SET
P-TKW1221	HOOK FOR BASKET	SET
P-TKW1209	KNEE PAD	EACH
P-TKW1212	KNOB WITH SCREW FOR HANDLE HEIGHT SEE SN REFERENCE BELOW	EACH
P-TKW1206	LEVER FOR FOLDING MECHANISM	EACH
P-TKW1223	PLATFORM - NEW POST	EACH
P-TKW1210	PLATFORM - OLD POST SEE SN REFERENCE BELOW	EACH
P-TKW1218	PLUG FOR MAIN FRAME AT THE END	EACH
P-TKW1402	POST LOWER STEERING	EACH
P-TKW1203	SCREW & WASHER FOR BRAKE	EACH
P40103	SCREW & NUT TO ATTACH BRAKE HANDLE	PAIR
P-TKW1217	SCREW FOR BRAKE ADJUSTMENT	PAIR

	Item Number TKW-12	
	THE KNEE WALKER	
	SERIAL # BEGINNING WITH "NT"	
ITEM NUMBER	DESCRIPTION	UOM
P-TKW1407	SCREW FOR LOWER STEERING POST	EACH
P-TKW1211	SCREW SET (8 PCS) FOR SEAT	SET
P-TKW1208	SCREW W/ KNOB FOR HANDLE BAR AND SEAT HT. ADJ	EACH
TKC-1041	SPACER FOR WHEEL BEARING	EACH
P-TKW1216	SPRING (PAIR) FOR BRAKE	PAIR
P-TKW1226	STEERING POST	EACH
P-TKW1219	STEERING POST SEE SN REFERENCE BELOW	EACH
P-TKW1204	WHEEL	SET
P-TKW1220	WRENCH (2 PCS)	SET

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>PRE NOV. 2016 MODELS</i>
TKW-12	SN BEFORE: NTTKW-12FK0251

SERIAL NUMBER (SN) REFERENCE	
ITEM NUMBER	<i>POST MAR. 2020 MODELS</i>
TKW-12	SN BEFORE: NTTKW-12KB0001


CARE • QUALITY • INNOVATION



THE NOVA PROMISE

CONTACT US!
We're here to help

 info@novajoy.com

 1 800 557 6682

novajoy.com