



Fall 2021

nova

Bathroom Safety



# Do You Know the BATHROOM HOT SPOTS?



We use the bathroom every day ... 365 days of each year, yet this room is often the most dangerous place in your home. Being proactive with effective product solutions is the key to fall prevention.


For less than a dollar a day, you can equip your bathroom with products to keep you safe. Make your bathroom happy and safe every day.



GET SOME LEVERAGE

## HOT SPOT #1


If you have trouble getting on and off a chair, then you are also having trouble getting on and off a toilet. An elevated toilet seat and hand rails can make all the difference.



GETTING IN & OUT OF THE TUB

## HOT SPOT #2

Getting in and out of the tub or shower can be tricky, and holding onto the sink or towel rack isn't safe. Using grab bars is a much better solution.



BATHING INDEPENDENTLY

## HOT SPOT #3

Our tub or shower is our sanctuary, but without the proper safety equipment, this is another bathroom hot spot. Bath seats, hand-held shower heads and bath mats are a must!

## VISIT US ON YouTube

NOVA's **YouTube** channel is jam packed with **great information, inspiration** and that **NOVAjoy**. We have **over 140 videos and counting** to educate, inform and inspire our customers.

Please do **Subscribe to our Channel** and feel free to give us a "thumbs up" or feedback. And, thank you for watching and sharing the joy!



Scan here to see our YouTube Channel!



## GRAB BARS

### Wall Grab Bars 8012-R to 8032-R



Chrome Ivory

- 1.5" wall stand-off
- 2.75" fixed flange
- 1" diameter tubing
- Knurled chrome grip (Chrome only)

ITEM #	DESCRIPTION
8012CH-R	12" Wall Grab Bar - Chrome
8016CH-R	16" Wall Grab Bar - Chrome
8018CH-R	18" Wall Grab Bar - Chrome
8024CH-R	24" Wall Grab Bar - Chrome
8032CH-R	32" Wall Grab Bar - Chrome
8012PC-R	12" Wall Grab Bar - Ivory
8016PC-R	16" Wall Grab Bar - Ivory
8018PC-R	18" Wall Grab Bar - Ivory
8024PC-R	24" Wall Grab Bar - Ivory
8032PC-R	32" Wall Grab Bar - Ivory

### Home Grab Bar 8609W-R

- Sleek, modern design with low profile
- Non-slip grip on bar promotes stability
- 250 lb weight capacity
- Only requires one wood stud to install and can be used throughout the home
- Best used in home danger zones such as bathrooms, doorways & stairwells
- Dimensions: 12.5"w x 2"d x 2.5"h
- Weight Capacity: 250 lb



### Luxury Wall Grab Bar 8618N-R

- 18" Length
- Nickel
- Corrosion resistant
- Textured surface for improved grip
- Hardware included
- 1.5" wall stand-off
- 3.25" fixed flange
- 1.25" diameter tubing



Textured surface for improved grip

## Suction Cup Balance Assist Bar


8212-R

- Installs quickly; perfect for home or travel use
- Provides added support when entering and leaving a bath tub or shower
- Easy to attach and remove
- Secure suction indicator displays green for secure and red for not secure
- No tools required
- Dimensions: 5" w x 3.5" d x 12" h




## 16" Suction Cup Balance Assist Bar

8216B-R / 8216N-R



8216B-R  
Bronze

8216N-R  
Nickel



- Green and red indicators help ensure proper installation
- Waved tubing promotes secure hand grip
- Strong, easy-to-use suction cup mechanism
- Use for balance assist only, not for body weight leverage
- Each suction cup requires only 5" x 5" of smooth, non-porous tile space
- Dimensions: 4.25" w x 3.25" d x 20.5" h

## Hand Held Shower Set

9301-R (2 Settings)

- Extra long 84" reinforced hose allows user to shower while seated
- On/Off switch allows for water flow control
- Includes wall mounting holder



9301-R

## Hand Held Shower Sets

9311-R & 9312-R (5 Settings)

- Extra long 84" reinforced hose allows user to shower while seated
- On/Off switch allows for water flow control
- Includes wall mounting holder
- 5 Settings: Regular / Shower Spray / Massage / Water Drip / Shower Stop



9311-R  
(White)

9312-R  
(Chrome)

## Tub Grab Bar

8220-R



- Provides added support when entering and leaving a bath tub or shower
- Made of 1" powder coated steel tubing
- Rubber contact points
- Clamp adjusts from 3" - 7" w
- Dimensions: 6" w x 11" d x 19" h
- Weight Capacity: 250 lb

## Holder for Hand Held Shower

9306-R



- Attaches to most bath seat arms (3/4" diameter tube) or legs (1" diameter tube)
- Securely holds hand held shower in place for safety and convenience

## Suction Cup Showerhead Holder

8213-R



- Showerhead is in easy reach when seated or bathing children
- Suction cup can be quickly removed and relocated as often as you like
- Convenient when bathing pets or cleaning out the tub
- Works on surfaces that are flat, smooth and non-porous such as tile, glass, mirror, fiberglass and acrylic
- Dimensions: 2.5" w x 4.25" d x 3.75" h

## Suction Cup Towel Hook

8214N-R (Nickel)



- Easily add shower accessories to your bathroom
- Attaches via powerful suction mechanism without any tools
- Stylish finish for a modern touch to any bath or shower
- Easy to move or reposition

## Bath Seat Basket

9605-R



- Holds bath essentials conveniently within reach to avoid bending, turning & twisting
- Easily attaches to the seat frame with included hardware
- Fits all standard size bath benches with 7/8" or 1" tube diameter
- Coated steel wire is rust resistant
- Dimensions: 12" w x 4.5" d x 4.25" h

## Bath Seat Cushion

9601W-R



- Adds comfort while bathing
- Non-slip Velcro strap prevents cushion from sliding
- Fits all standard sized bath seats
- Color: White
- Dimensions: 14.75" w x 8.5" d x 1" h

## Rubber Bath Mats


9350-R, 9351-R & 9351W-R

9350-R



- Attaches securely to tub surface and helps prevent slipping
- Anti-fungal treated
- 15.75" w x 27.5" l

9351-R



Also available in white (9351W-R)

## Shower Diverter Valve

9305-R



- Designed to divert water between your shower head and hand-held shower

**Bath Seat**

**9039 with Arms & U-Shaped Cutout**

- U-shaped cutout for improved hygiene
- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>16" - 21" | • <b>Leg Footprint:</b><br>High - 19.5"w x 18.5"d |
| • <b>Seat Dimensions:</b><br>15.75"w x 13"d   | • <b>Weight Capacity:</b> 300 lb                  |

**Heavy Duty Bath Seat**

**9051 with U-Shaped Cutout**

- U-shaped cutout for improved hygiene
- Adjustable height aluminum frame
- Latex-free rubber tips



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>15.75" - 21.75" | • <b>Leg Footprint:</b><br>High - 20"w x 18.5"d |
| • <b>Seat Dimensions:</b><br>19"w x 12.5"d          | • <b>Weight Capacity:</b> 500 lb                |

**Heavy Duty Bath Seat**

**9061 with Back & U-Shaped Cutout**

- U-shaped cutout for improved hygiene
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>15.75" - 21.75" | • <b>Leg Footprint:</b><br>High - 20"w x 18.5"d |
| • <b>Seat Dimensions:</b><br>19"w x 12.5"d          | • <b>Weight Capacity:</b> 500 lb                |

**Bath & Shower Seats**

**9010 & 9020**

- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9020 only)




- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>16" - 22" | • <b>Leg Footprint:</b><br>High - 20"w x 17.75"d |
| • <b>Seat Dimensions:</b><br>20"w x 12.25"d   | • <b>Weight Capacity:</b> 300 lb                 |

**Bath & Shower Seats**

**9050 & 9060 with U-Shaped Cutout**

- U-shaped design for cleaning purposes
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9060 only)



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>14.25" - 21.25" | • <b>Leg Footprint:</b><br>High - 20"w x 18"d |
| • <b>Seat Dimensions:</b><br>19"w x 12.75"d         | • <b>Weight Capacity:</b> 300 lb              |

**Bath & Shower Seats**

**9016 & 9026**

- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips
- Removable back (9026 only)



- |   |   |
|---|---|
| • <b>Seat Height:</b> 16" - 20.75"          | • <b>Leg Footprint:</b><br>High - 19.25"w x 18.25"d |
| • <b>Seat Dimensions:</b><br>15.75"w x 13"d | • <b>Weight Capacity:</b> 300 lb                    |
| • <b>Width Between Arms:</b> 18.75"         |   |



## Deluxe Bath & Shower Seats

### 9120-R & 9100-R

- Safety handle for extra support and balance can be attached to either side.
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips



Red Safety Handle: Item # 9201-R

- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>15" - 21" | • <b>Leg Footprint:</b><br>High - 17" w x 18.5" d<br>Low - 16.75" w x 17.25" d |
| • <b>Seat Dimensions:</b><br>19.5" w x 16" d  | • <b>Weight Capacity:</b> 300 lb   |

## Heavy Duty Bath & Shower Seats

### 9013 & 9023

- Adjustable height aluminum frame
- Built-in hand held shower holder
- Latex-free rubber tips
- Removable back (9023 only)



9013 9023

- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>17" - 23" | • <b>Leg Footprint:</b><br>High - 19.75" w x 18" d<br>Low - 18.5" w x 16.75" d |
| • <b>Seat Dimensions:</b><br>20" w x 12" d    | • <b>Weight Capacity:</b> 500 lb   |

## Bath Seat

### 9131-R with Back

- Tool-free assembly in under 5 minutes
- Ergonomic, oversized seat supports up to 500 lb
- Stylish and modern colors give a classic feel to your bathroom
- Slip-resistant rubber tips add stability
- Non-corrosive aluminum legs will not rust



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>16.5" - 20.5" | • <b>Leg Footprint:</b><br>High - 19.25" w x 16.5" d<br>Low - 18.5" w x 16" d |
| • <b>Seat Dimensions:</b><br>17" w x 14.75" d     | • <b>Weight Capacity:</b> 500 lb  |

## Bath & Shower Seats

### 9121-R & 9101-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Built-in hand held shower holder
- Adjustable height aluminum frame
- Latex-free rubber tips



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>17" - 22" | • <b>Leg Footprint:</b><br>High - 20" w x 17.5" d<br>Low - 18" w x 16.75" d |
| • <b>Seat Dimensions:</b><br>20" w x 12" d    | • <b>Weight Capacity:</b> 300 lb  |

## Bath & Shower Seats

### 9035-R & 9036-R

- Great for travel and easy to store
- Tool-free assembly in minutes
- Arms provide support for sitting & rising
- Adjustable height aluminum frame
- Latex-free rubber tips



- |   |  |
|---|--|
| • <b>Seat Height:</b> 16" - 21"               | • <b>Leg Footprint:</b><br>High - 19.75" w x 18" d<br>Low - 19" w x 16.75" d |
| • <b>Seat Dimensions:</b><br>15.75" w x 13" d | • <b>Weight Capacity:</b> 300 lb   |
| • <b>Width Between Arms:</b> 19"              |  |

## Bath Seat

### 9037 with Arms & U-Shaped Cutout

- Great for travel and easy to store
- U-shaped cutout for improved hygiene
- Tool-free assembly in minutes
- Arms provide support for sitting & standing
- Adjustable height aluminum frame
- Latex-free rubber tips



- |   |  |
|---|--|
| • <b>Seat Height:</b> 16" - 21"               | • <b>Leg Footprint:</b><br>High - 19.5" w x 18.5" d<br>Low - 18.75" w x 17.25" d |
| • <b>Seat Dimensions:</b><br>15.75" w x 13" d | • <b>Weight Capacity:</b> 300 lb   |
| • <b>Width Between Arms:</b> 18.75"           |  |

**Bath Stool**  
9006-R



- Great for travel and easy to store
- Tool-free assembly in minutes
- Adjustable height aluminum frame
- Latex-free rubber tips

• <b>Seat Height Adjustment:</b> 15.5" - 20.5"	• <b>Leg Footprint:</b> High - 15"w x 14.5"d
• <b>Seat Dimensions:</b> 12.5"w x 12.75"d	• <b>Weight Capacity:</b> 300 lb

**Bath Stool**  
9005



- Great for use in small areas
- Adjustable height aluminum frame
- Latex-free rubber tips

• <b>Seat Height Adjustment:</b> 13.75" - 20"	• <b>Leg Footprint:</b> High - 15.5"w x 14.75"d
• <b>Seat Dimensions:</b> 12.75"w x 12.75"d	• <b>Weight Capacity:</b> 250 lb

**Tub Shower Board**  
9200

- Width adjusts to attach securely to most bathtubs
- Built-in soap holder
- Rubber covered clamps resist slipping and protect surfaces from scratching
- Safety handle can be attached to either side for support



• <b>Seat Dimensions:</b> 31"w x 14.5"d	• <b>Weight Capacity:</b> 250 lb
• <b>Width Adjustment:</b> 17.75" - 28.5"	

**Shower Chair**  
9401 with U-Shaped Cutout



- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips

• <b>Seat Height Adjustment:</b> 18" - 20.75"	• <b>Leg Footprint:</b> High - 21.25"w x 17.25"d
• <b>Seat Dimensions:</b> 16"w x 15.75"d	• <b>Weight Capacity:</b> 250 lb

**Foldable Shower Chair**  
9400 with U-Shaped Cutout



- Great for travel/easy to store
- Folds and unfolds quickly
- U-shaped design for cleaning purposes
- Arms provide support for sitting and rising
- Adjustable height aluminum frame
- Latex-free rubber tips

• <b>Seat Height Adjustment:</b> 18" - 20.75"	• <b>Leg Footprint:</b> High - 21"w x 19.25"d
• <b>Seat Dimensions:</b> 16"w x 16"d	• <b>Weight Capacity:</b> 250 lb

**Wall Mounted Shower Seat**  
9404

- Mounts directly to shower wall
- Saves space by folding up when not in use
- Includes wall mount anchors
- Adjustable height aluminum frame
- Latex-free rubber tips



• <b>Seat Height Adjustment:</b> 19.75" - 22.25"	• <b>Overall Footprint:</b> High - 15.75"w x 13.75"d
• <b>Seat Dimensions:</b> 16"w x 16"d	• <b>Weight Capacity:</b> 250 lb

**Deluxe Transfer Bench****9070**

- Built-in soap holder and hand held shower holder
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>18.75" - 23.25" | • <b>Leg Footprint:</b><br>High - 33"w x 17.5"d |
| • <b>Seat Dimensions:</b><br>31.5"w x 14.5"d        | • <b>Weight Capacity:</b> 300 lb                |

**Standard Transfer Bench****9071**

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>17.25" - 21.25" | • <b>Leg Footprint:</b><br>High - 30"w x 20.75"d |
| • <b>Seat Dimensions:</b><br>27"w x 16.25"d         | • <b>Weight Capacity:</b> 300 lb                 |

**Heavy Duty Transfer Bench****9075**

- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>16.75" - 21.75" | • <b>Leg Footprint:</b><br>High - 29.5"w x 20.75"d |
| • <b>Seat Dimensions:</b><br>27"w - 16.25"d         | • <b>Weight Capacity:</b> 500 lb                   |

**Padded Transfer Bench****9080**

- Padded seat and back add comfort when bathing
- Converts easily for left or right hand transfers
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- |   |   |
|---|---|
| • <b>Seat Height Adjustment:</b><br>17.75" - 21.75" | • <b>Leg Footprint:</b><br>High - 29.75" x 20.5"d |
| • <b>Seat Dimensions:</b><br>27"w x 16"d            | • <b>Weight Capacity:</b> 300 lb                  |

**Portable Transfer Bench****9072-R**

- Great for travel and easy to store
- Tool-free assembly in minutes
- Converts easily for left or right hand transfers
- Latex-free rubber tips & suction cups



- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>19.25" - 23.25" | • <b>Leg Footprint:</b><br>High - 31.75"w x 22.5"d |
| • <b>Seat Dimensions:</b><br>26.5"w x 16.25"d       | • <b>Weight Capacity:</b> 300 lb                   |

**Padded Transfer Bench/Commode****9073**

- Converts easily for left or right hand transfers
- Removable seat for easy cleaning
- Bucket and lid included
- Adjustable height aluminum frame
- Latex-free rubber tips & suction tips



- |   |  |
|---|--|
| • <b>Seat Height Adjustment:</b><br>19" - 23.5" | • <b>Leg Footprint:</b><br>High - 36.5"w x 17.75"d |
| • <b>Seat Dimensions:</b><br>26.75"w x 14.5"d   | • <b>Weight Capacity:</b> 300 lb                   |



## Bedside Commode

8450



- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8.5" w x 9.5" d

- |   |   |
|---|---|
| • <b>Seat Height:</b> 17.5" - 23.5"           | • <b>Seat Dimensions:</b> 13.25" w x 14.25" d |
| • <b>Overall Dimensions:</b> 24" w x 19.75" d | • <b>Weight:</b> 11.5 lb                      |
| • <b>Width Between Arms:</b> 17.75"           | • <b>Weight Capacity:</b> 350 lb              |

## Folding Commode

8700-R



- Can be used beside the bed and over the toilet as a raised toilet seat with arms
- Ideal for travel: folds and unfolds in seconds
- Adjustable height steel frame
- Bucket, lid, and splash guard included
- Toilet seat opening: 8" w x 10" d

- |  |   |
|--|---|
| • <b>Seat Height:</b> 16.25" - 22.25"          | • <b>Seat Dimensions:</b> 13.75" w x 16.25" d |
| • <b>Overall Dimensions:</b> 22.5" w x 21.5" d | • <b>Weight:</b> 10 lb                        |
| • <b>Width Between Arms:</b> 18.5"             | • <b>Weight Capacity:</b> 300 lb              |

## Drop-Arm Commode

8900W



- Arms drop for easy patient transfer
- Bucket, lid, and splash guard included
- Available in Gray and Chrome
- Adjustable height steel frame
- Toilet seat opening: 7.75" w x 9.75" d

- |  |   |
|--|---|
| • <b>Seat Height:</b> 18.75" - 22.25"                    | • <b>Seat Dimensions:</b> 13.5" w x 16" d |
| • <b>Overall Dimensions:</b> 22.5" w x 20.75" - 22.75" d | • <b>Weight:</b> 19.6 lb                  |
| • <b>Width Between Arms:</b> 18"                         | • <b>Weight Capacity:</b> 300 lb          |

## Padded Drop-Arm Commode

8901W



- Arms drop for easy patient transfer
- Padded seat and backrest
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet seat opening: 8.25" w x 11" d

- |   |   |
|---|---|
| • <b>Seat Height:</b> 19.25" - 23.75"                 | • <b>Seat Dimensions:</b> 16.5" w x 16.5" d |
| • <b>Overall Dimensions:</b> 22.5" w x 21.25" - 23" d | • <b>Weight:</b> 20 lb                      |
| • <b>Width Between Arms:</b> 17.25"                   | • <b>Weight Capacity:</b> 300 lb            |

## Heavy Duty Commode

8580



- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet seat opening: 8" w x 10" d

- |  |   |
|--|---|
| • <b>Seat Height:</b> 19.75" - 23.5"             | • <b>Seat Dimensions:</b> 23" w x 18.5" d |
| • <b>Overall Dimensions:</b> 32.25" w x 23.25" d | • <b>Weight:</b> 22.25 lb                 |
| • <b>Width Between Arms:</b> 28.5"               | • <b>Weight Capacity:</b> 450 lb          |

## Heavy Duty Commode

8582



- Extra-large seat provides added support
- Bucket, lid, and splash guard included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25" w x 10.25" d

- |   |   |
|---|---|
| • <b>Seat Height:</b> 17.5" - 24.25"          | • <b>Seat Dimensions:</b> 23.25" w x 18.75" d |
| • <b>Overall Dimensions:</b> 27" w x 23.75" d | • <b>Weight:</b> 16.6 lb                      |
| • <b>Width Between Arms:</b> 23.75"           | • <b>Weight Capacity:</b> 500 lb              |

## Heavy Duty Drop-Arm Commode

8583

- Arms drop for easy patient transfer
- Extra-large seat provides added support
- Bucket and lid included
- Adjustable height steel frame
- Toilet Seat Opening: 8.25" w x 10.25" d



- |                                     |  |
|-------------------------------------|--|
| • Seat Height: 19.75" - 23.75"      | • Seat Dimensions: 23.25" w x 18.75" d |
| • Overall Dimensions: 30" w x 22" d | • Weight: 18.1 lb                      |
| • Width Between Arms: 25.75"        | • Weight Capacity: 500 lb              |

## Drop-Arm Transport Chair Commode

8805

- Arms drop for easy patient transfer
- 4" swivel wheels with locking casters
- Padded seat fits over commode and converts into transport chair
- Bucket & lid included
- Removable, swing away footrests
- Bucket can be removed without lifting user or raising seat
- Toilet Seat Opening: 8.75" w x 10" d
- NOT FOR USE IN SHOWER



- |  |                                       |
|--|---------------------------------------|
| • Seat Height: 18.75"                          | • Seat Dimensions: 17.75" w x 17.5" d |
| • Overall Dimensions: 21" w x 22.25" d x 38" h | • Weight: 22 lb                       |
| • Width Between Arms: 17.5"                    | • Weight Capacity: 250 lb             |

## Shower Chair/Commode

8800

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Toilet seat opening: 8.5" w x 10" d



- |  |                                       |
|--|---------------------------------------|
| • Seat Height: 22.25"                            | • Seat Dimensions: 15.5" w x 16.75" d |
| • Overall Dimensions: 22.75" w x 22.5" d x 34" h | • Weight: 13.25 lb                    |
| • Width Between Arms: 19.5"                      | • Weight Capacity: 250 lb             |

## Deluxe Shower Chair/Commode

8801

- Waterfall armrests and padded seat
- 4" swivel wheels with locking casters
- Bucket and lid included
- Removable, swing away footrests
- Toilet seat opening: 8.5" w x 10" d



- |  |                                  |
|--|----------------------------------|
| • Seat Height: 22.25"                          | • Seat Dimensions: 16" w x 16" d |
| • Overall Dimensions: 22" w x 24" d x 36.25" h | • Weight: 18.1 lb                |
| • Width Between Arms: 18.75"                   | • Weight Capacity: 250 lb        |

## Commode Liner

PA-1080

- Medical grade bag with super absorbent pad that quickly transitions body fluids into a gel, eliminating odors, spills & splashes
- Provides convenience & safety for Patient, Caregiver & Family
- No more hand-cleaning of soiled pail
- No handling of feces or urine
- Convenient drawstring closure for easier cleanup
- Pad absorbs 17 oz. of body fluids



## Padded Toilet Seat Risers

2629-R (2") & 2630-R (4")

- High density foam minimizes pressure on sensitive areas
- Attaches easily to toilet seat with four adjustable Velcro straps to reduce shifting
- Easy to clean vinyl surface
- Fits most toilet seats
- Dimensions:
  - 2629-R: 14" w x 15.75" d x 2" h
  - 2630-R: 14.5" w x 15.75" d x 4" h
- Weight Capacity: 250 lb



2629-R



2630-R

## Toilet Safety Rails

8200-R



- One piece mounting flange for easy width adjustment
- Rails bolt securely to the toilet
- Waterfall armrest

• **Arm Height Adjustment:**  
26.5" - 30.5"

• **Width Between Arms:**  
21" - 24"

• **Weight:** 3 lb

• **Weight Capacity:** 250 lb

## Toilet Support Rails

8201-R



- Provides added stability when sitting and rising
- Frame bolts directly to toilet
- Cushioned, easy to clean foam armrests
- Adjusts to fit standard and elongated toilets
- Easy to install

• **Arm Height:**  
11.5"

• **Width Between Arms:**  
17.25" - 19.75"

• **Overall Dimensions:**  
23.5"w x 16.25"d x 11.5"h

• **Weight Capacity:** 300 lb

## Toilet Safety Support Frame

8205-R



- Freestanding frame
- Slip-resistant base
- Padded foam handles
- For use with both standard and elongated toilets
- Includes magazine rack
- Easy, tool-free assembly

• **Overall Dimensions:**  
29.75"w x 19"d x 25.5"h

• **Handle Height:** 25.5"

• **Width Between Arms:** 20.5"

• **Weight Capacity:** 300 lb

## Toilet Seat Risers

8341-R & 8342-R



- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.75" height to the toilet seat
- 8341-R (Elongated)
- 8342-R (Standard)

• **Overall Dimensions:**  
14"w x 19.25"d x 3.75"h (8341-R)

• 13.75"w x 17.25"d x 3.75"h (8342-R)

• **Weight:** 2.4 lb (8341) & 2.3 lb (8342)

• **Aperture Dimensions:**  
10.25"w x 13"d (8341-R)

• 9.5"w x 11"d (8342-R)

• **Weight Capacity:** 300 lb

## Toilet Seat Risers

8343-R & 8344-R with Arms



- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Rubber pads help prevent seat from moving
- Adds 3.75" height to the toilet seat
- Padded arms
- 8343-R (Elongated)
- 8344-R (Standard)

• **Overall Dimensions:**  
21.25"w x 19.25"d x 9.75"h (8343-R)

• 21"w x 17.25"d x 9.75"h (8344-R)

• **Weight:** 3.7 lb (8343) & 3.55 lb (8344)

• **Aperture Dimensions:**  
10.25"w x 13"d (8343-R)

• 9.5"w x 11.25"d (8344-R)

• **Weight Capacity:** 300 lb

## Toilet Seat Risers (Hinged)

8345-R & 8346-R



- Hinged design makes cleaning easier and faster than before
- Allows the original seat and lid to be used
- Bolts directly to the toilet for added safety
- Adds 3.75" height to the toilet seat
- 8346-R (Elongated)
- 8345-R (Standard)

• **Overall Dimensions:**  
13.75"w x 17"d x 3.75"h (8345-R)

• 14.25"w x 19.25"d x 3.75"h (8346-R)

• **Weight:** 2.7 lb (8345) & 3.05 lb (8346)

• **Aperture Dimensions:**  
9.5"w x 11.25"d (8345-R)

• 10.25"w x 14"d (8346-R)

• **Weight Capacity:** 300 lb

## 5" Raised Toilet Seat

8340-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Non-locking
- Extra large seating surface



- Overall Dimensions: 17.5"w x 17.25"d x 6.25"h
- Aperture Dimensions: 7.5"w x 9.5"d
- Weight: 3 lb
- Weight Capacity: 300 lb

## 5" Locking Raised Toilet Seat

8350-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation



- Overall Dimensions: 16.5"w x 17.25"d x 7.25"h
- Aperture Dimensions: 7.75"w x 9"d
- Weight: 4.15 lb
- Weight Capacity: 300 lb

## 5" Locking Raised Toilet Seat

8351-R with Adjustable Width Detachable Arms

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Arms provide support for sitting and rising
- Front locking mechanism for easy installation
- Padded arms



- Overall Dimensions: 22.25"w x 17"d x 14.25"h
- Aperture Dimensions: 7.25"w x 8.5"d
- Width Between Arms: 17.75" - 20"
- Weight: 5.25 lb
- Weight Capacity: 300 lb

## 5" Locking Raised Toilet Seat

8352-R

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Front locking mechanism for easy installation



- Overall Dimensions: 20.5"w x 17.75"d x 7"h
- Aperture Dimensions: 8.25"w x 9.5"d
- Weight: 4.6 lb
- Weight Capacity: 300 lb

## 5" Locking Raised Toilet Seat

8353-R with Detachable Arms

- Adds 5" to the height of the toilet seat to aid in sitting and rising
- Arms provide support for sitting and rising
- Front locking mechanism for easy installation
- Padded arms



- Overall Dimensions: 20.5"w x 17.75"d x 12.75"h
- Aperture Dimensions: 8"w x 9.5"d
- Width Between Arms: 18"
- Weight: 5.75 lb
- Weight Capacity: 300 lb

## Easy Air Adjustable Toilet Seat Riser

2628A-R (2.25")

- Adjustable air foam combination
- Air release valve easily adjusts to your individual comfort
- Automatically re-inflates with no pump required
- Easily attaches to your toilet seat with loop straps
- Dimensions: 14"w x 15.5"d x 2.25"h
- Weight Capacity: 250 lb





# PREVENT FALLS

This is one epidemic we can cure



## IT'S NOT AN AGE THING

Fall related injuries in the bathroom are as common among those 18-25 as they are among those over 65. Let's make the bathroom safe for everyone!

Consult your physician before  
use of any medical products



as unique as you are

1-800-557-6682

[novajoy.com](http://novajoy.com)

Authorized NOVA Dealer