STAR HD

Parts Breakdown

For models: 4259BL 4259PL, 4259RD

Serial # Starting "HF"

HEAVY DUTY HD

400 lb Weight Capacity



CARE • QUALITY • INNOVATION

THE NOVA PROMISE

Revised on May 9, 2025

CONTACT US!
We're here to help

info@novajoy.com

1 800 557 6682

novajoy.com

An important note from

Many parts listed and/or pictured on this breakdown may appear to be identical to the part(s) required but may actually be different, due to slight dimensional changes.

Please verify that the serial number on your rolling walker corresponds with the serial number series denoted in the below parts list. The serial number is located on a sticker that is attached to the frame under the seat (see Figure 1).


If you need further assistance, please call our customer service team at (800) 557-6682.

Thank you!




Figure 1


Hand Brake Mechanism




Upper Assembly Components




Seat, Backrest & Bag



Lower Assembly Components



Fork, Wheel, & Brake Shoe Assembly



	Item Number 4259	
	STAR HD	
	SERIAL # BEGINNING WITH "HF"	
ITEM NUMBER	DESCRIPTION	UOM
HF-425901	AXLE & NUT FRONT	EACH
HF-425902	AXLE & NUT REAR	EACH
HF-425910	BACKREST	EACH
HF-425909	BAG	EACH
HF-425903	BEARING & SPACER WHEEL	SET
VT-33207	BEARING FOR SWIVEL FRONT FORK	EACH
P-CABLE	BRAKE CABLE (PAIR)	PAIR
HF-425907	BRAKE CABLE HOLDER	PAIR
HF-425809	BRAKE SHOE (SET)	SET
HF-420506	CLAMP AND SCREW (SET) FOR SEAT	SET
JN-423624	COVER FOR CABLE	EACH
P40013	CRIMP FOR STEEL CABLE (PAIR)	PAIR
HF-425904	FOLDING MECHANISM	EACH
JN-423615	GRIPS (PAIR)	PAIR
HF-425807	HAND BRAKE ASSEMBLY	EACH
HF-425912	HANDLE LEFT	EACH
HF-425913	HANDLE RIGHT	EACH
HF-425914	HANDLE/LEG ASSEMBLY LEFT	EACH
HF-425915	HANDLE/LEG ASSEMBLY RIGHT	EACH
HF-425804	HINGE BACKREST LEFT	EACH
HF-425820	HINGE BACKREST RIGHT	EACH
HF-425905	LEG FRONT LEFT ASSEMBLY	EACH
HF-425918	LEG FRONT RIGHT ASSEMBLY	EACH
HF-425906	LEG REAR LEFT	EACH
HF-425919	LEG REAR RIGHT	EACH
JN-420439	REFLECTOR (PAIR) FOR HAND BRAKE	PAIR
HF-425805	SCREW A-FRAME AND BACKREST HINGE	EACH

	Item Number 4259	
	STAR HD	
	SERIAL # BEGINNING WITH "HF"	
ITEM NUMBER	DESCRIPTION	UOM
C40006	SCREW TO ATTACH HAND BRAKE	EACH
HF-425908	SEAT PADDED	EACH
HF-425911	SILENCER FOR FRAME	EACH
JN-423911	STEM SCREW FOR FRONT FORK	EACH
JN-423626	TENSION SCREW	EACH
HF-425916	TOE KICK (PAIR) FOR LEG	PAIR
HF-425917	WHEEL 8"	EACH
HF-425820	BACKREST PADDED FOR 4259	EACH
HF-425921	WHEEL 8" FOR 4259 WB	EACH

STARTING STRING FOR HIGHLIGHTED PARTS:
HF4259SD0001**

CONTACT US!
We're here to help
info@novajoy.com
1 800 557 6682
novajoy.com

CARE • QUALITY • INNOVATION

THE NOVA PROMISE